Sartorius MA150
The compact class featuring maximum performance with minimum space requirements

For routine operation
A rugged design with low space requirements and easy operation are the major features of the MA150. Fully automatic drying of a sample until a constant weight is reached eliminates the need for programming an endpoint shutoff parameter. Twenty drying routines can be saved to give you the flexibility you need when the moisture content of additional, "out-of-the-ordinary" samples of material has to be measured.

Customizable and fast
Sartorius offers you a choice of two different moisture analyzers that cover diverse requirements on moisture measurements. Whichever heat source you opt for, both analyzers deliver results within just minutes. For temperature-sensitive samples, a ceramic heating element ensures especially gentle heating over the entire surface. The other choice, a CQR quartzglass heater, optimizes the analysis time even further, which is already ultrafast for the analyzer featuring the ceramic heater.

Application-specific solutions
Practical accessories round off the entire line-up of Sartorius moisture analyzers. These include, for instance, an in-use dust cover that is included with the standard equipment supplied and a special version of the moisture analyzer without openly accessible glass components in compliance with the stringent FDA and HACCP requirements that ban the use of glass in production.





