PCR FAST

TECHNICAL DATASHEET



MADE IN ITALY

TECHNICAL SPECIFICATIONS

PRODUCT SAFETY CABINET

PCR FAST | EN 61010-1, ISO 14644-1 with HEPA-Filter

Description	Unit	
General dimensions W×H×D	mm	727 × 960 × 620
Usable dimensions W×H×D	mm	700 × 550 × 600
Work opening	mm	200
Max. front opening	mm	500
Weight	kg	55
UV radiation flow rate	mW	170
Average service life of the UV lamp	h	8.000
UV light timer setting	h	4 (+/- 1 min accuracy)
UV light peak wavelength	nm	275
Noise level[1]	db(A)	<54
Light intensity	lux	>2100
Internal protection	A	4
Electrical connection		1Ph+E - 230V 50Hz
Power consumption[2]	kW	0,11 kW
Equipment		- HEPA filter + motor fan kit - Underframe for seated workstation - Additional power outlet

[1] Equipped with a HEPA filter kit.
[2] Clean filters, lighting activated, internal output load excluded.

DASITGROUP FAST ER

EN/Tec_PCR_FAST_EN/2026-05/Ed. 01 In this brochure all pictures and specifications are purely indicative and may vary without notice.